

Anticon 100 Wipes

100% no-run interlock knit polyester with knife cut edges

Anticon 100 wipes are made from 100% continuous filament polyester with a double-knit interlock construction. The clean-cut edges exhibit low particle and fiber counts to minimize contamination. Using patented textile technology, Anticon 100 achieves superior capacity with a quick sorption rate. Ideal for spill control, general use as well as applying and removing disinfectants and cleaning solutions.

Anticon 100 wipes are available in StandardWeight, HeavyWeight, and Sterile 100 fabrics. The Anticon Sterile 100 is made of the Anticon 100 Heavyweight material and gamma irradiated to a Sterility Assurance Level of 10⁻⁶ SAL which eliminates the introduction of biological contamination. It is ideal for critical pharmaceutical environments.



Features	Benefits
Two weight options	• Lighter basis weight provides superior value for general cleaning • Heavy basis weight achieves superior capacity with a quick sorption rate
100% continuous filament polyester	• Specially processed for superior cleanliness and sorbency
Cut edge	• Low fibres reduces risk of contamination
Gamma irradiated to a Sterility Assurance Level of 10 ⁻⁶ SAL	• Minimize the risk of contamination

Part No.	Description	Size	Packaging
495352-803	Anticon 100 StandardWeight Wipes Flat stacked	230 x 230 mm	75 per bag / 16 bags per case
495352-805	Anticon 100 StandardWeight Wipes Flat stacked	305 x 305 mm	50 per bag / 8 bags per case
492222-813	Anticon 100 HeavyWeight Wipes Flat stacked	230 x 230 mm	75 per bag / 16 bags per case
492222-951	Anticon 100 HeavyWeight Wipes Bulk	230 x 230 mm	75 per bag / 16 bags per case
492222-815	Anticon 100 HeavyWeight Wipes Flat stacked	305 x 305 mm	50 per bag / 8 bags per case
492222-808	Anticon 100 HeavyWeight Wipes Flat stacked	305 x 305 mm	25 per bag / 24 bags per case
492225-762	Anticon 100 HeavyWeight Wipes Bulk	305 x 305 mm	63 per bag / 5 bags per case
492226-939	Anticon Sterile 100 HeavyWeight Wipes Flat stacked	230 x 230 mm	20 per bag / 40 bags per case

Product Information	
Material	100% polyester
Construction	Interlock knit
Environment	ISO 3-8 Grade A/B for sterile, C/D for nonsterile



Technical Data			
	StandardWeight	HeavyWeight	
Attribute (units)	Typical Value	Typical Value	Test Method
Basis weight, nominal; (g/m ²)	120	140	Contec Method
Sorbent capacity; (mL/m ²)	339	492	IEST-RP-CC004.3, Sec. 8.1
Sorptive rate; (seconds)	<1	<1	
Non-volatile residue, NVR			IEST-RP-CC004.3, Sec. 7.1.1
In deionized water; (g/m ²)	ND	0.009	
In isopropyl alcohol; (g/m ²)	ND	0.003	
Specific ions			IEST-RP-CC004.3, Sec. 7.2.2
Sodium; (ppm)	0.184	0.022	
Chloride; (ppm)	0.030	0.007	
Particles, readily releasable			
Particles ≥ 0.5µm; (x10 ⁶ /m ²)	3.7	4.1	IEST-RP-CC004.3, Sec. 6.1.3
Fibres ≥ 100µm; (x 10 ³ /m ²)	1.413	1.70	IEST-RP-CC004.3, Sec. 6.2.2

Notes

- The data shown are typical values and should not be used as product specifications.
- Valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.
- Current and/or comparison data may be available. Please contact a Contec sales representative for details.

Test Methods

- CTM = Contec Test Method
- IEST-RP-CC004.3 = Evaluating Wiping Materials Used in Cleanroom and Other Controlled Environments, Institute of environmental Sciences and Technology, Rolling Meadows IL.

Packaging

Packaging Materials				
Outer bags	Low density polyethylene (LDPE) 			
Case	Corrugated fibreboard (PAP) 			
Packaging Configuration	EA/OB1	OB1/OB2	OB2/CS	EA/CS
495352-803	75	2	8	1,200
495352-805	50	2	4	400
492222-813	75	2	8	1,200
492222-951	75	2	8	1,200
492222-815	50	2	4	400
492222-808	25	3	8	600
492225-762	63	1	5	315
492226-939	20	5	8	800

EA = Each; OB1 = Outer Bag 1; OB2 = Outer Bag 2; CS = Case

Recycle Key

PET	
LDPE	
PAP	

For more information or to request a sample,
email infoeu@contecinc.com or phone +33 (0) 2 97 43 76 98.